

Anti-Prostaglandin E synthase 2 PTGES2 Antibody

Catalog Number: A06706

Overview

Product Name	Anti-Prostaglandin E synthase 2 PTGES2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Prostaglandin E synthase 2 PTGES2 Antibody catalog # A06706. Tested in ELISA,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8BWM0

Technical Details

Immunogen	KLH conjugated synthetic peptide derived from mouse PTGES2:281-384/384
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	ELISA: 1:5000-10000 IHC: 1:100-500 IF: 1:100-500

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Prostaglandin E synthase 2 PTGES2 Antibody

For Research Use Only. Not for use in diagnostic procedures.