

Anti-SPRR3 Antibody Picoband® HRP Conjugated

Catalog Number: A06729-2-HRP

About SPRR3

Predicted to enable structural molecule activity. Predicted to be involved in keratinocyte differentiation and wound healing. Located in Golgi apparatus and perinuclear region of cytoplasm.

Overview

Product Name	Anti-SPRR3 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9UBC9

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human SPRR3.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.