

## Anti-SUMO4 Antibody

Catalog Number: A06740


### Overview

Product Name	Anti-SUMO4 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SUMO4 Antibody catalog # A06740. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6EEV6


### Technical Details

Immunogen	Recombinant full length Human SUMO4.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

## Anti-SUMO4 Antibody (A06740) Images



Western blot analysis of SUMO4 polyclonal antibody



Immunohistochemistry (IHC) analysis of SUMO4 polyclonal antibody

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-SUMO4 Antibody

For Research Use Only. Not for use in diagnostic procedures.