

Anti-SUMO4 Antibody

Catalog Number: A06740-1

Overview

Product Name	Anti-SUMO4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SUMO4 Antibody catalog # A06740-1. Tested in WB,ICC/IF,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6EEV6

Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 FC: 1:50-1:100 ICC: 1:100-1:500

Anti-SUMO4 Antibody (A06740-1) Images




Figure 2. Immunocytochemistry staining of SUMO4 using Anti-SUMO4 Antibody (A06740-1).

ICC staining SUMO4 in 293 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde




Figure 1. Western blotting validation for Anti-SUMO4 Antibody A06740-1

Western blot analysis of SUMO4 on different lysates using anti-SUMO4 antibody at 1/1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SUMO4 Antibody