

Anti-HSL (Phospho-S853) LIPE Antibody

Catalog Number: A06762S853


Overview

Product Name	Anti-HSL (Phospho-S853) LIPE Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HSL (Phospho-S853) LIPE Antibody catalog # A06762S853. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q05469

Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser853 of human Hormone sensitive lipase.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000

Anti-HSL (Phospho-S853) LIPE Antibody (A06762S853) Images



Western blot analysis of phospho -Hormone sensitive lipase (S853) on human skeletal muscle lysates using anti-phospho -Hormone sensitive lipase (S853) antibody at 1/1,000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSL (Phospho-S853) LIPE Antibody

For Research Use Only. Not for use in diagnostic procedures.