

Anti-MARCKS-related protein MARCKSL1 Antibody

Catalog Number: A06798

Overview

Product Name	Anti-MARCKS-related protein MARCKSL1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MARCKS-related protein MARCKSL1 Antibody catalog # A06798. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49006

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human MARCKSL1
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-MARCKS-related protein MARCKSL1 Antibody (A06798) Images




Figure 1. Western blotting validation for Anti-MARCKS-related protein MARCKSL1 Antibody A06798

Western blot (WB) analysis of MARCKSL1 polyclonal antibody at 1:500 dilution

Lane1:HEK293T whole cell lysate

Lane2:Raw264.7 whole cell lysate

Lane3:PC12 whole cell lysate

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MARCKS-related protein MARCKSL1 Antibody