

## Anti-Aprataxin and PNK-like factor APLF Antibody

Catalog Number: A06828-2


### Overview

Product Name	Anti-Aprataxin and PNK-like factor APLF Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Aprataxin and PNK-like factor APLF Antibody catalog # A06828-2. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8IW19


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human APLF. AA range:82-131
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:40000 IF 1:50-200


## Anti-Aprataxin and PNK-like factor APLF Antibody (A06828-2) Images



Western Blot analysis of HT-29 cells using APLF Polyclonal Antibody



Western blot analysis of lysates from HT-29 cells, using APLF Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using APLF Antibody. The picture on the right is blocked with the synthesized peptide.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Aprataxin and PNK-like factor APLF Antibody

For Research Use Only. Not for use in diagnostic procedures.