

## Anti-DP2/TFDP2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A06850-3-Fluoro647

### About TFDP2

Transcription factor Dp-2 is a protein that in humans is encoded by the TFDP2 gene. The gene is a member of the transcription factor DP family. The encoded protein forms heterodimers with the E2F transcription factors resulting in transcriptional activation of cell cycle regulated genes. Alternative splicing results in multiple transcript variants.

### Overview

Product Name	Anti-DP2/TFDP2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q14188

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human DP2/TFDP2, identical to the related mouse sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DP2/TFDP2 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.