

## Anti-Integrin alpha-9/ITGA9 Antibody Picoband® Biotin Conjugated

Catalog Number: A06851-1-Biotin

### About ITGA9

Integrin alpha-9 is a protein that in humans is encoded by the ITGA9 gene. This gene encodes an alpha integrin. Integrins are heterodimeric integral membrane glycoproteins composed of an alpha chain and a beta chain that mediate cell-cell and cell-matrix adhesion. The protein encoded by this gene, when bound to the beta 1 chain, forms an integrin that is a receptor for VCAM1, cytotactin and osteopontin. Expression of this gene has been found to be upregulated in small cell lung cancers.

### Overview

Product Name	Anti-Integrin alpha-9/ITGA9 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q13797

### Technical Details

Immunogen	E.coli-derived human Integrin alpha-9/ITGA9 recombinant protein (Position: Q526-Q1035). Human ITGA9 shares 89.6% amino acid (aa) sequence identity with mouse ITGA9.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

---

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin alpha-9/ITGA9 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.