

Anti-NFR kappaB NFRKB Antibody

Catalog Number: A06867

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-NFR kappaB NFRKB Antibody
Reactive Species	Human, Monkey, Mouse
Description	Boster Bio Anti-NFR kappaB NFRKB Antibody catalog # A06867. Tested in WB, ELISA, IHC applications. This antibody reacts with Human, Mouse, Monkey.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6P4R8


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human NFRKB. AA range:463-512
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified


Suggested Dilutions

WB 1:500-2000
IHC 1:50-300
ELISA 1:2000-20000

Anti-NFR kappaB NFRKB Antibody (A06867) Images



Western Blot (WB) analysis of COS7 cells using NFRkappaB Polyclonal antibody.



Western Blot (WB) analysis of specific cells using NFRkappaB Polyclonal antibody.

[Anti-NFR kappaB NFRKB Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.