

Anti-MIEF1 Antibody Picoband® PE Conjugated

Catalog Number: A06887-1-PE

About MIEF1

Mitochondrial dynamic protein MID51 (MID51) also known as mitochondrial elongation factor 1 (MIEF1) or Smith-Magenis syndrome chromosome region candidate gene 7 protein-like (SMCR7L) is a protein that in humans is encoded by the SMCR7L gene. The SMCR7L gene codes for a protein that has been called MiD51/MIEF1 and shown to regulate mitochondrial fission by interacting with the proteins Drp1 and FIS1.

Overview

Product Name	Anti-MIEF1 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IP, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NQG6

Technical Details

Immunogen	E.coli-derived human MIEF1 recombinant protein (Position: S189-T463).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MIEF1 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.