

## Anti-Casein Kinase I Gamma 1 (H3) CSNK1G1 Antibody

Catalog Number: A06904-1


### Overview

Product Name	Anti-Casein Kinase I Gamma 1 (H3) CSNK1G1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Casein Kinase I Gamma 1 (H3) CSNK1G1 Antibody catalog # A06904-1. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9HCP0


### Technical Details

Immunogen	Synthetic peptide, corresponding to the N-terminal of Human Casein Kinase Igamma1.
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200


## Anti-Casein Kinase I Gamma 1 (H3) CSNK1G1 Antibody (A06904-1) Images



Western blot (WB) analysis of Casein Kinase I gamma 1 (H3) pAb at 1:500 dilution Lane1:3T3-L1 whole cell lysate(40ug) Lane2:C6 whole cell lysate(40ug) Lane3:HEK293T whole cell lysate(40ug) Lane4:MCF-7 whole cell lysate(40ug) Lane5:A549 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of Casein Kinase I gamma 1 (H3) pAb in paraffin-embedded human brain tissue.



Western blot (WB) analysis of Casein Kinase I gamma 1 (H3) pAb at 1:500 dilution Lane1:MCF-7 whole cell lysate(40ug) Lane2:A549 whole cell lysate(40ug) Lane3:Hela whole cell lysate(40ug) Lane4:EC9706 whole cell lysate(40ug) Lane5:A2780 whole cell lysate(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Casein Kinase I Gamma 1 (H3) CSNK1G1 Antibody

For Research Use Only. Not for use in diagnostic procedures.