

Anti-ENTPD5 Antibody Picoband® HRP Conjugated

Catalog Number: A06908-2-HRP

About ENTPD5

The protein encoded by this gene is similar to E-type nucleotidases (NTPases)/ecto-ATPase/apyrases. NTPases, such as CD39, mediate catabolism of extracellular nucleotides. ENTPD5 contains 4 apyrase-conserved regions which is characteristic of NTPases.

Overview

| | |
|----------------------|---|
| Product Name | Anti-ENTPD5 Antibody Picoband® HRP Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | WB, IHC, ELISA |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | O75356 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human ENTPD5 recombinant protein (Position: D95-H428). Human ENTPD5 shares 90.1% and 88.9% amino acid (aa) sequence identity with mouse and rat ENTPD5, respectively. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ENTPD5 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.