

Anti-CD297 ART4 Antibody

Catalog Number: A06913


Overview

Product Name	Anti-CD297 ART4 Antibody
Reactive Species	Human
Description	Boster Bio Anti-CD297 ART4 Antibody catalog # A06913. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q93070


Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 181-230
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	IHC-P, 1:50-300

Anti-CD297 ART4 Antibody (A06913) Images



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD297 ART4 Antibody