

Anti-USP39 Antibody

Catalog Number: A06922

Overview

| | |
|----------------------|--|
| Product Name | Anti-USP39 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-USP39 Antibody catalog # A06922. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q53GS9 |

Technical Details

| | |
|---------------------|---|
| Form | Liquid |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | <p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB 1:1000 - 1:2000</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p> |

Anti-USP39 Antibody (A06922) Images




Figure 1. Western blotting validation for Anti-USP39 Antibody A06922

Western blot analysis of extracts of various cells

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-USP39 Antibody