

Anti-ARHGEF9 Antibody

Catalog Number: A06931-4

About ARHGEF9

ARHGEF9, also named as KIAA0424, acts as guanine nucleotide exchange factor (GEF) for CDC42. ARHGEF9 promotes formation of GPHN clusters. This antibody is specific to ARHGEF9.

Overview

Product Name	Anti-ARHGEF9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ARHGEF9 Antibody catalog # A06931-4. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O43307

Technical Details

Immunogen	Peptide
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ARHGEF9 Antibody

For Research Use Only. Not for use in diagnostic procedures.