

Anti-STX1B Antibody (N-term)

Catalog Number: A06932

About STX1B

Potentially involved in docking of synaptic vesicles at presynaptic active zones. May mediate Ca²⁺-regulation of exocytosis acrosomal reaction in sperm (By similarity).


Overview

Product Name	Anti-STX1B Antibody (N-term)
Reactive Species	Human
Description	Boster Bio Anti-STX1B Antibody (N-term) (Catalog # A06932). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P61266

Technical Details

Immunogen	This STX1B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human STX1B.
Predicted Reactive Species	Bovine, Mouse, Rat, Sheep
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

Anti-STX1B Antibody (N-term) (A06932) Images



STX1B Antibody (N-term) western blot analysis in SH-SY5Y cell line lysates (35ug/lane). This demonstrates the STX1B antibody detected the STX1B protein (arrow).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-STX1B Antibody (N-term)

For Research Use Only. Not for use in diagnostic procedures.