

Anti-TAUT SLC6A6 Antibody

Catalog Number: A06936

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-TAUT SLC6A6 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TAUT SLC6A6 Antibody catalog # A06936. Tested in ELISA, IHC-P, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P31641


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human SLC6A6. AA range:561-610
Isotype	Rabbit IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using a epitope-specific immunogen.

Suggested Dilutions

WB 1:500-1:2000
IHC 1:100-1:300
ELISA 1:40000
IF 1:50-200

Anti-TAUT SLC6A6 Antibody (A06936) Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using SLC6A6 Antibody. The picture on the right is blocked with the synthesized peptide.

[Anti-TAUT SLC6A6 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.