

Anti-MCTS1 Antibody

Catalog Number: A06981-2

About MCTS1

MCTS1 is an anti-oncogene that plays a role in cell cycle regulation. It can reduce cell doubling time and anchor-dependent growth. The G1 transition time and the duration of the G1 / S transition are shortened. It can also promote the release of deacylated tRNA and mRNA from recycled 40S subunits following ABCE1-mediated dissociation of post-termination ribosomal complexes into subunits (PMID: 10440924).

Overview

Product Name	Anti-MCTS1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MCTS1 Antibody catalog # A06981-2. Tested in WB, IHC, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9ULC4

Technical Details

Immunogen	E.coli-derived human MCTS1 recombinant protein (Position: K3-K178).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MCTS1 Antibody

For Research Use Only. Not for use in diagnostic procedures.