

Anti-CYP2C8/9/18/19 (M136) CYP2C18 Antibody

Catalog Number: A07005-1


Overview

Product Name	Anti-CYP2C8/9/18/19 (M136) CYP2C18 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CYP2C8/9/18/19 (M136) CYP2C18 Antibody catalog # A07005-1. Tested in IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P33260/P10632/P11712/P33261

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 100-150 of Human CYP2C8.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200 IF: 1:50-1:200

Anti-CYP2C8/9/18/19 (M136) CYP2C18 Antibody (A07005-1) Images



Immunohistochemistry (IHC) analyzes of CYP2C8/9/18/19 (M136) pAb in paraffin-embedded liver cancer tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CYP2C8/9/18/19 (M136) CYP2C18 Antibody

For Research Use Only. Not for use in diagnostic procedures.