

Anti-Multimerin-1 MMRN1 Antibody

Catalog Number: A07012


Overview

Product Name	Anti-Multimerin-1 MMRN1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Multimerin-1 MMRN1 Antibody catalog # A07012. Tested in ELISA applications. This antibody reacts with Human,Mouse,Rat.
Application	ELISA
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13201


Technical Details

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1083-1228 of human MMRN1 (NP_031377.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	ELISA: 1:5000-1:10000

Anti-Multimerin-1 MMRN1 Antibody (A07012) Images



Western blot analysis of various lysates, using MMRN1 Rabbit pAb at 1:1000 dilution.



Immunofluorescence analysis of U2OS cells using MMRN1 antibody. Blue: DAPI for nuclear staining.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Multimerin-1 MMRN1 Antibody

For Research Use Only. Not for use in diagnostic procedures.