

Anti-TBL1X Monoclonal Antibody

Catalog Number: A07014


Overview

Product Name	Anti-TBL1X Monoclonal Antibody
Reactive Species	Human, Monkey
Description	Boster Bio Anti-TBL1X Monoclonal Antibody catalog # A07014. Tested in WB, ELISA applications. This antibody reacts with Human, Monkey.
Application	ELISA, WB
Clonality	Monoclonal 4G1
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O60907

Technical Details

Immunogen	Purified recombinant fragment of human TBL1X expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-TBL1X Monoclonal Antibody (A07014) Images



Western Blot analysis using TBL1X Monoclonal Antibody against HEK293 (1) cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TBL1X Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.