

## Anti-MC5R Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A07019-3-Fluoro550

### About MC5R

Predicted to enable hormone binding activity and melanocortin receptor activity. Predicted to be involved in adenylate cyclase-activating G protein-coupled receptor signaling pathway and regulation of metabolic process. Predicted to be active in cytoplasm and plasma membrane. Is expressed in skin. Orthologous to human MC5R (melanocortin 5 receptor).

### Overview

Product Name	Anti-MC5R Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P41149

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of mouse MC5R. Mouse MC5R shares 95.5% amino acid (aa) sequence identity with both human and rat MC5R.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-MC5R Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.