

Anti-RPL18 (H151) Antibody

Catalog Number: A07057

Overview

Product Name	Anti-RPL18 (H151) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-RPL18 (H151) Antibody catalog # A07057. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q07020

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 130-180 of Human Ribosomal Protein L18.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:10-1:100

