

Anti-Krueppel-like factor 13 KLF13 Antibody

Catalog Number: A07076-1


Overview

Product Name	Anti-Krueppel-like factor 13 KLF13 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Krueppel-like factor 13 KLF13 Antibody catalog # A07076-1. Tested in ELISA, IHC, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y2Y9


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human KLF13 around the non-acetylation site of Lys166. AA range:131-180
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-1:300 ELISA 1:20000 IF 1:50-200


Anti-Krueppel-like factor 13 KLF13 Antibody (A07076-1) Images




Western blotting validation for Anti-Krueppel-like factor 13 KLF13 Antibody A07076-1 Western Blot (WB) analysis of L929 cells using KLF13 polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunohistochemistry validation of KLF13 using Anti-Krueppel-like factor 13 KLF13 Antibody (A07076-1). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Pancreas



Immunohistochemistry validation of KLF13 using Anti-Krueppel-like factor 13 KLF13 Antibody (A07076-1). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Pancreas



Immunohistochemistry validation of KLF13 using Anti-Krueppel-like factor 13 KLF13 Antibody (A07076-1). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Lung

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.