

## Anti-SLC6A12 Antibody Picoband® Biotin Conjugated

Catalog Number: A07083-1-Biotin

### About SLC6A12

Enables gamma-aminobutyric acid:sodium:chloride symporter activity. Involved in gamma-aminobutyric acid transport and monocarboxylic acid transport. Predicted to be located in basolateral plasma membrane. Predicted to be active in cell projection and plasma membrane.

### Overview

Product Name	Anti-SLC6A12 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P48065

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human SLC6A12.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC6A12 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.