

## Anti-SLC6A12 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A07083-1-Fluoro550

### About SLC6A12

Enables gamma-aminobutyric acid:sodium:chloride symporter activity. Involved in gamma-aminobutyric acid transport and monocarboxylic acid transport. Predicted to be located in basolateral plasma membrane. Predicted to be active in cell projection and plasma membrane.

### Overview

Product Name	Anti-SLC6A12 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P48065

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human SLC6A12.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-SLC6A12 Antibody - Fluoro550**

**For Research Use Only. Not for use in diagnostic procedures.**