

Anti-Integrin Alpha E ITGAE Antibody

Catalog Number: A07091-1


Overview

Product Name	Anti-Integrin Alpha E ITGAE Antibody
Reactive Species	Human
Description	Boster Bio Anti-Integrin Alpha E ITGAE Antibody catalog # A07091-1. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P38570

Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	IHC: 1:50-1:200

Anti-Integrin Alpha E ITGAE Antibody (A07091-1) Images



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-Integrin alpha E antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin Alpha E ITGAE Antibody

For Research Use Only. Not for use in diagnostic procedures.