

Anti-EFTUD1/EFL1 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A07093-1-Fluoro594

About EFL1

Enables GTPase activity and ribosome binding activity. Involved in GTP metabolic process and cytosolic ribosome assembly. Predicted to be part of ribonucleoprotein complex. Predicted to be active in cytosol. Implicated in Shwachman-Diamond syndrome.

Overview

Product Name	Anti-EFTUD1/EFL1 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q7Z2Z2

Technical Details

Immunogen	E.coli-derived human EFTUD1/EFL1 recombinant protein (Position: E188-D932).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-EFTUD1/EFL1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.