

## Anti-EFTUD1/EFL1 Antibody Picoband® HRP Conjugated

Catalog Number: A07093-1-HRP

### About EFL1

Enables GTPase activity and ribosome binding activity. Involved in GTP metabolic process and cytosolic ribosome assembly. Predicted to be part of ribonucleoprotein complex. Predicted to be active in cytosol. Implicated in Shwachman-Diamond syndrome.

### Overview

Product Name	Anti-EFTUD1/EFL1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q7Z2Z2

### Technical Details

Immunogen	E.coli-derived human EFTUD1/EFL1 recombinant protein (Position: E188-D932).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.