

## Anti-Serum amyloid A-4 protein SAA4 Antibody

Catalog Number: A07115-1


### Overview

Product Name	Anti-Serum amyloid A-4 protein SAA4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Serum amyloid A-4 protein SAA4 Antibody catalog # A07115-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35542

### Technical Details

Immunogen	A synthetic peptide of human SAA4
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:2000-1:5000 IHC: 1:50-1:200

## Anti-Serum amyloid A-4 protein SAA4 Antibody (A07115-1) Images



Western blot analysis of extracts of various cell lines, using SAA4 antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Serum amyloid A-4 protein SAA4 Antibody

For Research Use Only. Not for use in diagnostic procedures.