

Anti-UBAP2L Antibody

Catalog Number: A07183-1

About UBAP2L

The UBA domain-containing protein ubiquitin associated protein 2-like (UBAP2L), also known as NICE4, is a ubiquitously expressed 1,087 amino acid protein. UBAP2L is preferentially expressed in hematopoietic stem cell (HSC) populations and interacts with BMI1, a well-known determinant of HSC function (1). It is essential for the activity of long-term repopulating HSCs but not the BMI1-dependent repression of the Ink4a/Arf locus (1).

Overview

Product Name	Anti-UBAP2L Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-UBAP2L Antibody (Catalog # A07183-1). Tested in ELISA, WB, ICC, IF applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, ICC, WB
Clonality	Polyclonal
Formulation	UBAP2L Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	UBAP2L antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Host	Rabbit
Uniprot ID	Q14157


Technical Details

Immunogen	UBAP2L antibody was raised against a 19 amino acid peptide from near the amino terminus of human UBAP2L. The immunogen is located within amino acids 200 - 250 of UBAP2L.
Predicted Reactive Species	Bovine
Cross Reactivity	UBAP2L antibody is human and mouse reactive. At least four isoforms of UBAP2L are known to exist; this antibody will detect all four isoforms. This antibody is predicted to not cross-react with UBAP2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	UBAP2L Antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	UBAP2L antibody can be used for the detection of UBAP2L by Western blot at 1 - 2 ug/mL. Antibody

can also be used for immunocytochemistry at 10 ug/ml.

Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.


Anti-UBAP2L Antibody (A07183-1) Images



Western blot analysis of UBAP2L in HeLa cell lysate with UBAP2L antibody at 1 ug/ml.



Immunocytochemistry of UBAP2L in HeLa cells with UBAP2L antibody at 10 ug/ml.



Immunofluorescence of UBAP2L in HeLa cells with UBAP2L antibody at 20 ug/ml. Red: UBAP2L Antibody (A07183-1)
Blue: DAPI staining

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBAP2L Antibody

For Research Use Only. Not for use in diagnostic procedures.