

Anti-Ras-related protein Rab-39B RAB39B Antibody

Catalog Number: A07186

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Ras-related protein Rab-39B RAB39B Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Ras-related protein Rab-39B RAB39B Antibody catalog # A07186. Tested in WB applications. This antibody reacts with Human,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96DA2

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human RAB39B
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Ras-related protein Rab-39B RAB39B Antibody (A07186) Images



Western blot (WB) analysis of RAB39B polyclonal antibody at 1:500 dilution Lane1:MCF-7 whole cell lysate Lane2:H9C2 whole cell lysate

Anti-Ras-related protein Rab-39B RAB39B Antibody

For Research Use Only. Not for use in diagnostic procedures.