

Anti-ACADL Antibody

Catalog Number: A07206


Overview

Product Name	Anti-ACADL Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ACADL Antibody catalog # A07206. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P28330


Technical Details

Immunogen	Recombinant fusion protein of human ACADL(NP_001599.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

Anti-ACADL Antibody (A07206) Images



Western blot analysis of ACADL expression in SHSY5Y (A), mouse liver (B), mouse heart (C), rat spinal cord (D) whole cell lysates.



Immunohistochemical analysis of ACADL staining in rat ovary formalin fixed paraffin embedded tissue section.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ACADL Antibody

For Research Use Only. Not for use in diagnostic procedures.