

## Anti-PKA I beta reg PRKAR1B Antibody

Catalog Number: A07227


### Overview

Product Name	Anti-PKA I beta reg PRKAR1B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PKA I beta reg PRKAR1B Antibody catalog # A07227. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P31321


### Technical Details

Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human FMNL3
Predicted Reactive Species	Canine, Equine, Guinea Pig, Rabbit, Zebrafish
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: WB, 1:500-1:2000

## Anti-PKA I beta reg PRKAR1B Antibody (A07227) Images



Western Blot (WB) analysis of COLO205 cells using PKA I beta reg Polyclonal antibody.



Western Blot (WB) analysis of specific cells using PKA I beta reg Polyclonal antibody.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PKA I beta reg PRKAR1B Antibody