

Anti-Ubiquitin-like protein 4A UBL4A Antibody

Catalog Number: A07314


Overview

Product Name	Anti-Ubiquitin-like protein 4A UBL4A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Ubiquitin-like protein 4A UBL4A Antibody catalog # A07314. Tested in WB,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P11441


Technical Details

Immunogen	Recombinant fusion protein of human UBL4A(NP_055050.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:200-1:1000 IP: 1:50-1:200

Anti-Ubiquitin-like protein 4A UBL4A Antibody (A07314) Images



Western blot analysis of extracts of various cell lines, using UBL4A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunoprecipitation analysis of 300ug extracts of HeLa cells using 3ug UBL4A antibody. Western blot was performed from the immunoprecipitate using UBL4A antibody at a dilution of 1:1000.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ubiquitin-like protein 4A UBL4A Antibody

For Research Use Only. Not for use in diagnostic procedures.