

## Anti-Mucin-20 MUC20 Antibody

Catalog Number: A07372-1


### Overview

Product Name	Anti-Mucin-20 MUC20 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Mucin-20 MUC20 Antibody catalog # A07372-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8N307

### Technical Details

Immunogen	Recombinant full length Human MUC20.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE)
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-Mucin-20 MUC20 Antibody (A07372-1) Images



Western blot (WB) analysis of MUC20 polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate(40ug) Lane2:The kidney tissue lysate of Mouse(40ug) Lane3:The kidney tissue lysate of Rat(40ug) Lane4:A549 whole cell lysate(40ug) Lane5:L02 whole cell lysate(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mucin-20 MUC20 Antibody

For Research Use Only. Not for use in diagnostic procedures.