

Anti-ADAM23 Antibody

Catalog Number: A07382


Overview

Product Name	Anti-ADAM23 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ADAM23 Antibody catalog # A07382. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75077


Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human ADAM23.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-ADAM23 Antibody (A07382) Images



Immunohistochemistry (IHC) analyzes of ADAM23 pAb in paraffin-embedded human rectum carcinoma tissue at 1:50, showing strong cytoplasmic staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.



Western blot (WB) analysis of ADAM23 polyclonal antibody at 1:500 dilution Lane1: MCF-7 whole cell lysate
Lane2: Raw264.7 whole cell lysate

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ADAM23 Antibody

For Research Use Only. Not for use in diagnostic procedures.