

Anti-METTL14 Antibody

Catalog Number: A07389-3

About METTL14

Methyltransferase like 14 is a protein that in humans is encoded by the METTL14 gene. METTL14 is part of a methyltransferase complex required for the common N6-methyladenosine (m6A) base modification that is enriched near stop codons and in long exons of mRNAs.

Overview

Product Name	Anti-METTL14 Antibody
Reactive Species	Human
Description	Boster Bio Anti-METTL14 Antibody catalog # A07389-3. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9HCE5

Technical Details

Immunogen	E.coli-derived human METTL14 recombinant protein (Position: Q111-R456).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-METTL14 Antibody

For Research Use Only. Not for use in diagnostic procedures.