

Anti-Cytokeratin 4 KRT4 Antibody

Catalog Number: A07410


Overview

Product Name	Anti-Cytokeratin 4 KRT4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cytokeratin 4 KRT4 Antibody catalog # A07410. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P19013

Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 200-260
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-2000 ELISA 1:10000-20000

Anti-Cytokeratin 4 KRT4 Antibody (A07410) Images



Western blot analysis of SKOV3 293T lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

1 Publications Citing This Product

1. PubMed ID: 22096478, Conjunctival Reconstruction with Progenitor Cell-Derived Autologous Epidermal Sheets in Rhesus Monkey

Visit bosterbio.com/anti-cytokeratin-4-antibody-a07410-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytokeratin 4 KRT4 Antibody

For Research Use Only. Not for use in diagnostic procedures.