

## Anti-TMX1 Antibody

Catalog Number: A07462-2

### About TMX1

TMX1, one of the few transmembrane members of the family, is the first vascular thiol isomerase with antithrombotic function. It forms functional complexes with the ER lectin calnexin and preferentially intervenes during maturation of cysteine-containing, membrane-associated proteins while ignoring the same cysteine-containing ectodomains if not anchored at the ER membrane. As such, TMX1 is the first example of a topology-specific client protein redox catalyst in living cells (PMID: 30655304). Covalent modification with AMS increases the molecular mass of originally oxidized species by ~0.5 kDa per thiol group, which allows enhanced separation from the reduced form in electrophoresis (PMID: 29123984). TMX1 has 2 bands produced by redox with the MW of 30-32 kDa.

### Overview

Product Name	Anti-TMX1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TMX1 Antibody catalog # A07462-2. Tested in WB, IHC, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9H3N1

### Technical Details

Immunogen	E.coli-derived human TMX1 recombinant protein (Position: D204-T277).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Immunohistochemistry, 1:50-400  
ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-TMX1 Antibody

For Research Use Only. Not for use in diagnostic procedures.