

Anti-Forkhead box protein I1 FOXI1 Antibody

Catalog Number: A07512-1


Overview

Product Name	Anti-Forkhead box protein I1 FOXI1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Forkhead box protein I1 FOXI1 Antibody (Catalog# A07512-1). Tested in IHC-P, IF (paraffin), WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC-P, WB, IF (paraffin)
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q12951

Technical Details

Immunogen	A peptide derived from human FOXI1. Immunogen sequence location: 191-240
Predicted Reactive Species	Bovine, Canine, Chicken, Primate, Sheep, Xenopus, Zebrafish
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>Western Blot, 1:500-2000</p> <p>Immunohistochemistry, 1:100 - 1:300</p> <p>ELISA, 1:5000</p> <p>Not yet tested in other applications.</p>

Anti-Forkhead box protein I1 FOXI1 Antibody (A07512-1) Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using FOXI1 Antibody. The picture on the right is blocked with the synthesized peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Forkhead box protein I1 FOXI1 Antibody