

Anti-TFIP11 Antibody (N-term)

Catalog Number: A07513-1

About TFIP11

TFIP11 is a nuclear speckle-localized protein that may play a role in spliceosome disassembly in Cajal bodies (Stanek et al., 2008 [PubMed 18367544]).


Overview

Product Name	Anti-TFIP11 Antibody (N-term)
Reactive Species	Human
Description	Boster Bio Anti-TFIP11 Antibody (N-term) (Catalog # A07513-1). Tested in IHC-P, WB application(s). This antibody reacts with Human.
Application	IHC-P, WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UBB9


Technical Details

Immunogen	This TFIP11 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 42-70 amino acids from the N-terminal region of human TFIP11.
Predicted Reactive Species	Monkey
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000 IHC-P: 1:25

Anti-TFIP11 Antibody (N-term) (A07513-1) Images



TFIP11 Antibody (N-term) western blot analysis in HeLa and K562 cell line and rat spleen tissue lysates (35ug/lane). This demonstrates the TFIP11 antibody detected the TFIP11 protein (arrow).



Immunohistochemical analysis of paraffin-embedded H. liver section using TFIP11 Antibody (N-term). A07513-1 was diluted at 1:25 dilution. A undiluted biotinylated goat polyvalent antibody was used as the secondary, followed by DAB staining.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TFIP11 Antibody (N-term)

For Research Use Only. Not for use in diagnostic procedures.