

## Anti-DPF2 (E174) Antibody

Catalog Number: A07556-2


### Overview

Product Name	Anti-DPF2 (E174) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DPF2 (E174) Antibody catalog # A07556-2. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92785

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 140-188 of Human DPF2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-DPF2 (E174) Antibody (A07556-2) Images



Western blot (WB) analysis of DPF2 (E174) polyclonal antibody at 1:1000 dilution Lane1:THP-1 whole cell lysate(37ug) Lane2:MCF-7 whole cell lysate(48ug) Lane3:RAW264.7 whole cell lysate(57ug) Lane4:NIH3T3 whole cell lysate(48ug) Lane5:H9C2 whole cell lysate(37ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-DPF2 (E174) Antibody

For Research Use Only. Not for use in diagnostic procedures.