

Anti-Collectin-11 COLEC11 Antibody

Catalog Number: A07583


Overview

Product Name	Anti-Collectin-11 COLEC11 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Collectin-11 COLEC11 Antibody catalog # A07583. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BWP8


Technical Details

Immunogen	Recombinant protein of human COLEC11.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:200-1:2000 IHC: 1:50-1:100


Anti-Collectin-11 COLEC11 Antibody (A07583) Images



Western blot analysis of extract of various cells, using COLEC11 antibody.



Immunohistochemistry of paraffin-embedded rat brain using COLEC11 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human prostate using COLEC11 antibody at dilution of 1:100 (40x lens).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Collectin-11 COLEC11 Antibody

For Research Use Only. Not for use in diagnostic procedures.