

Anti-SRRM2 Antibody

Catalog Number: A07592


Overview

Product Name	Anti-SRRM2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SRRM2 Antibody catalog # A07592. Tested in IHC, IF applications. This antibody reacts with Human, Mouse.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UQ35

Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 2000-2080
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	SRRM2 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:50-300 IF 1:50-200

Anti-SRRM2 Antibody (A07592) Images



Immunohistochemical analysis of paraffin-embedded human Small intestinal stromal tumor. 1, Tris-EDTA, pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200 (4° overnight). 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SRRM2 Antibody

For Research Use Only. Not for use in diagnostic procedures.