

## Anti-UDP-galactose translocator SLC35A2 Antibody

Catalog Number: A07601

### Overview


Product Name	Anti-UDP-galactose translocator SLC35A2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-UDP-galactose translocator SLC35A2 Antibody catalog # A07601. Tested in ICC/IF applications. This antibody reacts with Human.
Application	IF, ICC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P78381

### Technical Details

Immunogen	Recombinant fusion protein of human SLC35A2(NP_001027460.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	ICC/IF: 1:50-1:200

## Anti-UDP-galactose translocator SLC35A2 Antibody (A07601) Images

---



Immunofluorescence analysis of HeLa cells using SLC35A2 antibody . Blue: DAPI for nuclear staining.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UDP-galactose translocator SLC35A2 Antibody

For Research Use Only. Not for use in diagnostic procedures.