

## Anti-FNTB Antibody

Catalog Number: A07620


### Overview

Product Name	Anti-FNTB Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-FNTB Antibody (Catalog# A07620). Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49356

### Technical Details

Immunogen	A peptide derived from human FNTB. Immunogen sequence location: 255 - 304
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	WB 1:500 - 1:2000. ELISA: 1:40000. Not yet tested in other applications.

## Anti-FNTB Antibody (A07620) Images



Western blot analysis of lysates from RAW264.7 and NIH/3T3 cells, using FNTB Antibody. The lane on the right is blocked with the synthesized peptide.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-FNTB Antibody

For Research Use Only. Not for use in diagnostic procedures.