

Anti-NMT2 Antibody

Catalog Number: A07621


Overview

Product Name	Anti-NMT2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-NMT2 Antibody catalog # A07621. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60551


Technical Details

Immunogen	Recombinant full length Human NMT2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 IF: 1:50-1:200


Anti-NMT2 Antibody (A07621) Images



Western blot analysis of extracts of various cell lines, using NMT2 antibody.



Western blot analysis of extracts of various cell lines, using NMT2 antibody.



Immunofluorescence analysis of U2OS cell using NMT2 antibody. Blue: DAPI for nuclear staining.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NMT2 Antibody

For Research Use Only. Not for use in diagnostic procedures.