

Anti-RGC32/RGCC Antibody

Catalog Number: A07644-2

About RGCC

RGC32, also named as C13orf15, modulates the activity of cell cycle-specific kinases. RGC32 enhances CDC2 activity. RGC32 may contribute to the regulation of the cell cycle. And it inhibits growth of glioma cells by promoting arrest of mitotic progression at the G2/M transition. The antibody can recognize these 2 isoforms of RGC32.

Overview

Product Name	Anti-RGC32/RGCC Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-RGC32/RGCC Antibody catalog # A07644-2. Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9H4X1

Technical Details

Immunogen	Peptide
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RGC32/RGCC Antibody

For Research Use Only. Not for use in diagnostic procedures.